

PROVIDING A PREMIUM  
QUALITY FULL RANGE OF

# Medtopia Syringes

## Luer Lock & Luer Slip



CODE	DESCRIPTION	SIZE	UOM
M0203-06	Medtopia Syringe, luer lock	1ml	100pcs/box, 32boxes/ctn
M0203-01	Medtopia Syringe, luer lock	3ml	100pcs/box, 30boxes/ctn
M0203-02	Medtopia Syringe, luer lock	5ml	100pcs/box, 24boxes/ctn
M0203-03	Medtopia Syringe, luer lock	10ml	100pcs/box, 16boxes/ctn
M0203-04	Medtopia Syringe, luer lock	20ml	50pcs/box, 16boxes/ctn
M0206-06	Medtopia Syringe, luer slip	1ml	100pcs/box, 32boxes/ctn
M0206-01	Medtopia Syringe, luer slip	3ml	100pcs/box, 30boxes/ctn
M0206-02	Medtopia Syringe, luer slip	5ml	100pcs/box, 24boxes/ctn
M0206-03	Medtopia Syringe, luer slip	10ml	100pcs/box, 16boxes/ctn
M0206-04	Medtopia Syringe, luer slip	20ml	50pcs/box, 16boxes/ctn
M0206-05	Medtopia Syringe, luer slip	60ml	100pcs/box, 16boxes/ctn

**+ PRODUCT FEATURES:**

- Compatible with any standard hypodermic needles
- Sterilization Method : Ethylene oxide (ETO)
- Silicone coated barrel for a smooth plunger action
- DEHP free & Latex free

**+ CLASSIFICATION:**  
MDA registration No: GA5922523-120179

**+ PRODUCT MATERIAL:**  
High transparent medical grade PP (barrel) & PE (plunger)


### Luer Lock

### Luer Slip



**STORAGE RECOMMENDATION:**  
Store in cool place between 10C - 39C. Opened boxes should be avoid exposing to sunlight of fluorescent lighting to prevent discoloration.

**SHELF LIFE:**  
Expiring 5 years from manufactured date.



**LIFE VALUE**

**DISCLAIMER:**

- Do not use if unit package is damaged.
- Dispose after use to avoid risk of infections.

